

Dilation & Achalasia Balloon



Dilation Balloon

Catheter Outer Ø (mm/fr)	Working Length	Balloon Length	Balloon Ø	Recommended Pressure (ATM)	Burst Pressure (ATM)	Box Quantity	Product Code
2.4 / 7	250cm	40cm	6mm	8	10	1	CJ-JHY-BD-06-40-250
2.4 / 7	250cm	40cm	8mm	8	10	1	CJ-JHY-BD-08-40-250
2.4 / 7	250cm	40cm	10mm	8	10	1	CJ-JHY-BD-10-40-250
2.4 / 7	250cm	60cm	12mm	5	6	1	CJ-JHY-BD-12-60-250
2.4 / 7	250cm	60cm	14mm	5	6	1	CJ-JHY-BD-14-60-250
2.4 / 7	250cm	60cm	16mm	5	6	1	CJ-JHY-BD-16-60-250
2.4 / 7	250cm	80cm	16mm	5	6	1	CJ-JHY-BD-16-80-250
2.4 / 7	250cm	80cm	18mm	4	5	1	CJ-JHY-BD-18-80-250
2.4 / 7	250cm	80cm	20mm	4	5	1	CJ-JHY-BD-20-80-250

Achalasia Balloon

Catheter Outer Ø (mm/fr)	Working Length	Balloon Length	Balloon Ø	Recommended Pressure (ATM)	Burst Pressure (ATM)	Box Quantity	Product Code
6.7 / 20	90cm	80cm	30	2	3	1	CJ-JHY-BD-30-80-90
6.7 / 20	90cm	80cm	35	2	3	1	CJ-JHY-BD-35-80-90
6.7 / 20	90cm	80cm	40	2	3	1	CJ-JHY-BD-40-80-90

