

GI UK Medical – Data sheet
Dilation Balloons
 3 Diameter



Compliant Dilation Balloons

Catheter Outer Ø	Working Length	Recommended Pressure	Balloon Diameter	Balloon Length	Box Quantity	Product Code
2.4mm	240cm	3 / 4.5 / 6	10 / 11 / 12	55cm	1	NLS/BD3-10-55-240-GW
2.4mm	240cm	3 / 4.5 / 6	12 / 13.5 / 15	55cm	1	NLS/BD3-12-55-240-GW
2.4mm	240cm	3 / 4.5 / 6	15 / 16.5 / 18	55cm	1	NLS/BD3-15-55-240-GW
2.4mm	240cm	3 / 4.5 / 6	18 / 19 / 20	55cm	1	NLS/BD3-18-55-240-GW
2.4mm	240cm	3 / 4.5 / 6	12 / 13.5 / 15	80cm	1	NLS/BD3-12-80-240-GW
2.4mm	240cm	3 / 4.5 / 6	15 / 16.5 / 18	80cm	1	NLS/BD3-15-80-240-GW
2.4mm	240cm	3 / 4.5 / 6	18 / 19 / 20	80cm	1	NLS/BD3-18-80-240-GW

