


# Injection Needles



- Ergonomic handle to control needle movements
- Transparent catheter for liquid flow visualisation

Supplier	Product Code	Needle		Catheter		QOM	FAL Code
		Diameter (G / mm)	Length (mm)	Diameter (mm)	Length (cm)		
MDD	IN-23-19-4-230	19 / 1.0	4	2.4	230		FAL23697
MDD	IN-23-19-5-230	19 / 1.0	5	2.4	230		
MDD	IN-23-19-6-230	19 / 1.0	6	2.4	230		FAL23656
MDD	IN-23-22-4-160	22 / 0.7	4	2.4	160		
MDD	IN-23-22-4-230	22 / 0.7	4	2.4	230		FAL23521
MDD	IN-23-22-5-160	22 / 0.7	5	2.4	160		
MDD	IN-23-22-5-230	22 / 0.7	5	2.4	230		
MDD	IN-23-22-6-230	22 / 0.7	6	2.4	230		FAL23517
MDD	IN-23-23-4-160	23 / 0.6	4	2.4	160		
MDD	IN-23-23-4-230	23 / 0.6	4	2.4	230		FAL23831
MDD	IN-23-23-5-160	23 / 0.6	5	2.4	160		
MDD	IN-23-23-5-230	23 / 0.6	5	2.4	230	10	FAL23679
MDD	IN-23-23-6-230	23 / 0.6	6	2.4	230		FAL23574
MDD	IN-23-25-4-230	25 / 0.5	4	2.4	230		FAL23898
MDD	IN-23-25-5-230	25 / 0.5	5	2.4	230		FAL23630
MDD	IN-23-25-6-230	25 / 0.5	6	2.4	230		FAL23572
Pentax	IN-18-120-2205-A	22 / 0.5	120	1.8	180		
Pentax	IN-23-180-2204-A	22 / 0.4	180	2.4	180		
Pentax	IN-23-230-1904-A	19 / 0.4	230	2.4	230		
Pentax	IN-23-230-2205-A	22 / 0.5	230	2.4	230		
Pentax	IN-23-230-2504-A	25 / 0.4	230	2.4	230		
Pentax	IN-23-230-2506-A	25 / 0.6	230	2.4	230		