

GI UK Medical – Data sheet
Monofilament Snares



Hybrid Hot and Cold Monofilament Snares

Snare Width	Wire Ø	Catheter Ø	Catheter Length	Box Quantity	Product Code
MONOFILAMENT HOT SNARES					
10mm	0.3mm	2.3mm	230cm	10	CJ-ADR-23-230-MFO10
15mm	0.3mm	2.3mm	230cm	10	CJ-ADR-23-230-MFO15
20mm	0.3mm	2.3mm	230cm	10	CJ-ADR-23-230-MFO20
MONOFILAMENT TWO-STAGE					
15 & 25mm	0.3mm	2.3mm	230cm	10	CJ-DADR-23-230-025B
MULTI-STAGE ROTATABLE HYBRID SNARES					
10mm	0.3mm	2.4mm	230cm	10	OVA-HYB-24-10-230
15mm	0.3mm	2.4mm	240cm	10	OVA-HYB-24-15-230

