

RETRIEVAL NETS

GI UK Medical offers an all-encompassing range of retrieval nets for endoscopic use

RETRIEVAL NETS

Using the right retrieval net gives confidence in the capture & removal of foreign bodies. Our range offers versatility to deal with varying indications.



PRODUCT INFO

Oval and octagonal shape options offer varying degrees of stiffness to suit the environment and pathology being retrieved.

- The Blue outer wire clearly highlights location to improve visibility and accuracy of net placement.
- The translucent net enhances visibility.
- Our Memory Nylon net allows stretch to gather multiple tissue samples.
- Rotatable options for ease of positioning in tortuous positions.

STANDARD OVAL NETS

- Net cell spacing, guarantees visibility and tensile strength, whilst preserving sample integrity by mitigating sample crush.



OCTAGON NETS

- Memory points on the octagonal net increases rigidity and maintains device shape delivery increased durability.

FOREIGN BODY RETRIEVAL NETS					
Product Code	Net Opening		Catheter		Rotatable
	Width (mm)	Length (mm)	Diameter (mm)	Length (cm)	
STANDARD OVAL NETS					
FBN-S30-25-230-SU	20	30	2.5	230	-
FBN-S50-25-230-SU	25	50	2.5	230	-
FBN-S60-25-230-SU	30	60	2.5	230	-
FBN-S75-25-230-SU	40	70	2.5	230	-
NLS/POLRN-26-240-25	25	47	2.6	240	-
NLS/POLRN-26-240-32	32	52	2.6	240	-
OCTAGON NETS					
FBN-O35-25-230-SU	35	70	2.5	230	-
FBN-OR35-25-230-SU	35	70	2.5	230	Yes



43 Park Place, Leeds. LS1 2RY

Email: info@giukmedical.com

Web: www.giukmedical.com

Tel: 0113 513 4870