



# POLYP TRAPS

Our range of polyp traps are designed to simplify retrieval and identification of polyp

# P-TRAP

Polyp traps are integral to collection of small lesion pathology. Losing or damaging tissue samples has a detrimental impact on patients & procedures. GI UK Medical have a range of polyp traps to suit all clinical requirements.

## SOFT CATCH TRAYS

The soft catch tray system permits effective retrieval of pathology without compromising suction and loss of procedure time. The perforated drawers reduce risk of biomaterial exposure. The snap-lock system ensures drawers are securely closed to the trap, reducing risk of drawer handling error.



## MAGNIFYING GLASS DESIGN

Large draw system & magnifying glass design, allows clear identification of polyp retrieval.

- Quick and easy set up to aide theatre process flow.
- Integrity of pathology maintained with soft catch drawer design

Product Code	Description	FAL Code
VTL/PT-02	2 Drawer Polyp Trap	FAL23511



## SINGLE CHAMBER POLYP TRAP

- Single chamber trap offers an economic alternative for cases where limited pathology is to be collected.
- The clear transparent bubble case permits unobstructed view to confirm sample has been collected.

Product Code	Description	FAL Code
VTL/PT-01	1 Chamber Polyp Trap	FAL23721

## 5 CHAMBER POLYP TRAP

- 5-chamber rotatable trap permits multiple samples to be collected without loss of suction.
- 5 clearly marked chambers to match polyp excision location.
- Transparent chambers enables identification of polyp.
- Intuitive design allows simple and fast set up to support theatre flow

Product Code	Description	FAL Code
PT-01-05	5 Chamber Polyp Trap	FAL23720
VTL/PT-05	5 Chamber Polyp Trap	FAL23660





43 Park Place, Leeds. LS1 2RY

Email: [info@giukmedical.com](mailto:info@giukmedical.com)

Web: [www.giukmedical.com](http://www.giukmedical.com)

Tel: 0113 513 4870