

# Samarit Rollbord® for Lateral patient transfer

The SAMARIT Rollbord® (Rollboard) is the best aid for healthcare professionals to transfer patients safely and without damaging their backs in all areas of care and patient care in hospitals.

The Samarit „Rollbord“ is highly efficient, requiring only 15% of the patient's body weight to either push or pull the patient from one lateral surface to the other. One staff member can easily move an average sized patient with minimal effort.

Discover our 3 types of Rollboards:



	HIGHTEC	PROFESSIONAL	EcoLite
High quality product	yes	yes	2/3
Heavy duty use	yes	yes	no
Lowest friction (100kg)	15%	18%	20%
Safety endings	yes	yes	no
Antimicrobial towel	yes	yes	yes
Flame retardant	yes	yes	yes
Antistatic materials	yes	yes	yes
Light weighth constr.	yes	yes	yes
MRI combatible	yes	yes	yes
Samarit no-lift system	yes	yes	yes
Long life materials	yes	yes	yes
High workmanship	yes	yes	yes
Foldable models	yes	yes	yes
Numbers of models	5	4	3
Longlivity (years)	10	8	5
Easy to clean (wipes)	yes	yes	yes
Disposable covers	yes	yes	yes
Warranty years	3	2	2
SWL (safe weight limit)	< 310	< 180	< 160

## Hightec Rollbord - Standard



Model-No.:	440.0100
UDI/GTIN:	7640144460013
Use:	operating theatre, radiology, intensive care and special wards as well as in general care
Size:	178 x 50 cm / 89 x 50 cm foldable
Weight:	2.7 kg
Weight Limit:	240 kg

## Hightec Rollbord - ICU-Bariatric



Model-No.:	440.0200
UDI/GTIN:	7640144460020
Use:	Intensive care, general care
Size:	178 x 50 cm
Weight:	3.8 kg
Weight Limit:	240 kg

## Hightec Rollbord - Mini-Gyn



Model-No.:	440.0300
UDI/GTIN:	7640144460037
Use:	urology, gynecology, operating theatre, wake-up stations, emergencies
Size:	91 x 50 cm
Weight:	1.4 kg
Weight Limit:	180 kg

## Hightec Rollbord - Surgibord



Model-No.:	440.0400
UDI/GTIN:	7640144460044
Use:	operating theatre, radiology, MRI, CT, emergencies, general care as well as in homes
Size:	120 x 40 cm (Not foldable)
Weight:	1.5 kg
Weight Limit:	180 kg

## Hightec Rollbord - MRI/CT/AMBU



Model-No.:	440.0500
UDI/GTIN:	7640144460051
Use:	radiology, ambulance, emergency to transfer of injured patients on / from stretcher
Size:	178 x 40 cm / 87 x 40 cm foldable
Weight:	2.5 kg
Weight Limit:	240 kg



## Professional Rollbord - Standard



Model-No.:	440.2100
UDI/GTIN:	7640144460143
Use:	operating theatre, special wards as well as in general care
Size:	180 x 50 cm / 89 x 50 cm foldable
Weight:	2.9 kg
Weight Limit:	180 kg

## Professional Rollbord - MRI



Model-No.:	440.2300
UDI/GTIN:	7640144460150
Use:	radiology, MRI / CT
Size:	180 x 42 cm / 89 x 42 cm foldable
Weight:	2.3 kg
Weight Limit:	180 kg

## Professional Rollbord - AMBULANCE



Model-No.:	440.2400
UDI/GTIN:	7640144460167
Use:	Operating theatre, emergencies, ambulance as well as in homes
Size:	165 x 42 cm / 81 x 42 cm foldable
Weight:	1.9 kg
Weight Limit:	150 kg



## Ecolite Rollbord - ECOLITE 90



Model-No.:	440.1090
UDI/GTIN:	7640144460341
Use:	urology, gynecology, in emergencies and general care.
Size:	90 x 50 cm
Weight:	1.3 kg
Weight Limit:	150 kg

## Ecolite Rollbord - ECOLITE 120



Model-No.:	440.1120
UDI/GTIN:	7640144460334
Use:	radiology, emergencies as well as in homes
Size:	130 x 40 cm / 65 x 40 cm foldable
Weight:	1.4 kg
Weight Limit:	160 kg

## Ecolite Rollbord - ECOLITE 180



Model-No.:	440.1180
UDI/GTIN:	7640144460327
Use:	radiology, emergencies, general wards as well as in homes
Size:	180 x 50 cm / 89 x 50 cm foldable
Weight:	2.8 kg
Weight Limit:	160 kg

## Accessories - Dispenser for PE-Covers



The Rollbord Cover Dispenser holds one roll and includes a mounting hook for one Rollbord. The disposable covers are a specific length and are perforated at the tear point.

Size: 15 x 66 x 14 cm / Wall mount

## Accessories - Disposable Cover



To protect Rollbords from blood and dirt.

Model-No.: SAM 440.1000

UDI/GTIN: 7640144460471

1 box contains 4 rolls of 50 pieces (100m), perforated