



**DUOMED**

**Duomed Swiss AG**

Grenzstrasse 5a, 6214 Schenkon

Tel. +41 (0)41 510 07 00

[www.duomed.com](http://www.duomed.com)

# Endomyocardial Biopsy

- **JAWZ™ Endomyocardial Biopsy (EMB) Forceps**
- **Endomyocardial Biopsy Accessories**



## JAWZ™ Endomyocardial Biopsy (EMB) Forceps

JAWZ offers a clean-cut approach to better endomyocardial biopsy samples. It is designed to provide higher quality specimens with fewer passes required. The portfolio consists of disposables that maintain sharp dual cutting heads and have great design flexibility.

### Dual Cutting Head Yields Better Samples



- Unique cutting-edge design allows superior surface alignment to provide a clean cut with clean margins
- Precision-machined from a proprietary alloy to remain sharp throughout the procedure

### Enhanced Physician Experience

- Designed to deliver higher quality samples with reduced crushed artifact
- Precision-molded to provide a solid "instrument" feel

### Greater Flexibility & Ease of Use

- Three-ringed handle fully rotates and maintains your desired thumb position, staying in place
- Available in Stainless Steel or PTFE Coated, and Straight, Pre-Curved, and Maxi-Curved configurations



## Box of 1

Catalogue #	Jaw Size	Jaw OD (French)	Jaw Volume	Length	Curve	Coating	Formability	Fenestration
190090	1.5mm	5F	1.80mm <sup>3</sup>	50cm	Straight	None	Formable	Fenestrated
190095	1.5mm	5F	1.80mm <sup>3</sup>	105cm	Straight	None	Formable	Fenestrated
190051	1.8mm	6F	2.46mm <sup>3</sup>	50cm	Maxi-Curved	FEP	Formable	Non-Fenestrated
190060	1.8mm	6F	2.46mm <sup>3</sup>	50cm	Straight	None	Formable	Non-Fenestrated
190065	1.8mm	6F	2.46mm <sup>3</sup>	50cm	Straight	None	Non-Formable	Non-Fenestrated
190066	1.8mm	6F	2.46mm <sup>3</sup>	70cm	Straight	FEP	Formable	Non-Fenestrated
190070	1.8mm	6F	2.46mm <sup>3</sup>	105cm	Straight	None	Formable	Non-Fenestrated
190075	1.8mm	6F	2.46mm <sup>3</sup>	105cm	Straight	None	Non-Formable	Non-Fenestrated
190080	2.2mm	7F	4.97mm <sup>3</sup>	50cm	Pre-Curved	FEP	Formable	Non-Fenestrated
190081	2.2mm	7F	4.97mm <sup>3</sup>	50cm	Maxi-Curved	FEP	Formable	Non-Fenestrated
190085	2.2mm	7F	4.97mm <sup>3</sup>	105cm	Straight	None	Formable	Non-Fenestrated
190086	2.2mm	7F	4.97mm <sup>3</sup>	105cm	Straight	FEP	Formable	Non-Fenestrated
190020	2.4mm	7.5F	5.43mm <sup>3</sup>	50cm	Straight	None	Formable	Non-Fenestrated
190030	2.4mm	7.5F	5.43mm <sup>3</sup>	50cm	Pre-Curved	FEP	Formable	Non-Fenestrated
190031	2.4mm	7.5F	5.43mm <sup>3</sup>	50cm	Maxi-Curved	FEP	Formable	Non-Fenestrated
190040	2.4mm	7.5F	5.43mm <sup>3</sup>	105cm	Straight	None	Formable	Non-Fenestrated

## Endomyocardial Biopsy Accessories

### Access™ Percutaneous Introducer Kit w/Double-ended Wire (Box of 10)

Catalogue #	Sheath			Guidewire	
	Size	Length	Hub Colour	Diameter	Length
496004	4F	11cm	Red	.035"	45cm
496005	5F	11cm	Grey	.038"	45cm
496006	6F	11cm	Green	.038"	45cm
496007	7F	11cm	Orange	.038"	45cm
496008	8F	11cm	Blue	.038"	45cm
496009	9F	11cm	Black	.038"	45cm

### OneStic Micro-Access Introducer Kits (Box of 10)

Catalogue #	Sheath			Guidewire			Needle		
	Size	Length	Hub Colour	Diameter	Length	Core/Tip Material	Gauge	Length	Echogenic
497830	5F	11cm	Grey	.018"	45cm	SS/SS	21G	7cm	No
497840	6F	11cm	Green	.018"	45cm	SS/SS	21G	7cm	No
497841	6F	11cm	Green	.018"	40cm	NT/PD	21G	7cm	No
497861	7F	11cm	Orange	.018"	45cm	SS/SS	21G	7cm	No
497862	8F	11cm	Blue	.018"	45cm	SS/SS	21G	7cm	No

SS = Stainless Steel NT = Nitinol PD = Palladium