



DUOMED

Duomed Swiss AG

Grenzstrasse 5a, 6214 Schenkon

Tel. +41 (0)41 510 07 00

www.duomed.com

Guidewires

- **Guidewires - WORKER™**
- **Guidewires - Fixed Core Straight**
- **Guidewires - Fixed Core J Tip**
- **Guidewires - Lunderquist**
- **Guidewires - POINTER™**
- **Guidewires - Specialty**



Guidewires - WORKER™

WORKER guidewires are manufactured with a precoating technique offering a safe and flexible guidewire. The PTFE coating is used to reduce friction. All guidewires are packed in a circular dispenser facilitating easy flushing of the guidewire. Low friction due to PTFE coating; malleable tip; colour-coded guidewire straightener holds the wire in the dispenser and provides easy introduction in the puncture needle or catheter.



Standard - Straight, w/Fixed Core (Box of 10)


Catalogue #	Diameter	Length	Description	Soft Tip
110135080	.035"	80cm	Standard Straight	3.5cm
110135150	.035"	150cm	Standard Straight	3.5cm
110135180	.035"	180cm	Standard Straight	3.5cm
110135300	.035"	300cm	Standard Straight	3.5cm
110135400	.035"	400cm	Standard Straight	3.5cm

Standard - w/ 3mm J Curve, Fixed Core (Box of 10)

Catalogue #	Diameter	Length	Description	Soft Tip
110535080	.035"	80cm	Standard J Curve	3.5cm
110535150	.035"	150cm	Standard J Curve	3.5cm
110535180	.035"	180cm	Standard J Curve	3.5cm
110535260	.035"	260cm	Standard J Curve	3.5cm

Worker Amplatz

The Worker Amplatz guidewire is the right combination of stiffness, durability, and rail support.



Amplatz - Straight (Box of 10)

Catalogue #	Diameter	Length	Description	Soft Tip
114135080	.035"	80cm	Amplatz Straight	3.5cm
114135150	.035"	150cm	Amplatz Straight	3.5cm
114135180	.035"	180cm	Amplatz Straight	3.5cm
114135260	.035"	260cm	Amplatz Straight	3.5cm
114035090	.035"	90cm	Amplatz Straight	7.5cm
114035150	.035"	150cm	Amplatz Straight	7.5cm
114035180	.035"	180cm	Amplatz Straight	7.5cm
114035260	.035"	260cm	Amplatz Straight	7.5cm

Amplatz - 3mm J Curve (Box of 10)

Catalogue #	Diameter	Length	Description	Soft Tip
114235090	.035"	90cm	Amplatz J Curve	7.5cm
114235150	.035"	150cm	Amplatz J Curve	7.5cm

Guidewires - Fixed Core J Tip

PTFE Coated Stainless Steel (Box of 10)

Catalogue #	Diameter	Length	Description	Tip Curve Radius
395420	.025"	150cm	Original; Our original wire with traditional feel and performance	1.5mm
395787	.035"	150cm	Original; Our original wire with traditional feel and performance	1.5mm
395523	.021"	150cm	Original; Our original wire with traditional feel and performance	3mm
395559	.025"	150cm	Original; Our original wire with traditional feel and performance	3mm
395600	.032"	150cm	Original; Our original wire with traditional feel and performance	3mm
388282	.035"	150cm	Classic; A stiffer body wire for increased pushability	3mm
390282	.035"	150cm	XTB; A stiffer body wire with a shorter taper for increased pushability	3mm
393182	.035"	150cm	Zebra Tip; Partially masked (from PTFE coating) J tip for increased tip extendability	3mm
393282	.035"	150cm	Flexible; Our softer tip wire for increased maneuverability	3mm
395282	.035"	150cm	Original; Our original wire with traditional feel and performance	3mm
395284	.038"	150cm	Original; Our original wire with traditional feel and performance	3mm

PTFE Coated Stainless Steel Double-ended (Box of 10)

Catalogue #	Diameter	Length	Tip Curve Radius
395919	.035"	150cm	Str & 1.5mm

Guidewires - Fixed Core Straight

Mullins Tip Deflector Wires - Stainless Steel (Box of 10)

Catalogue #	Diameter	Length
395945	.018"	100cm

Stainless Steel (Box of 10)

Catalogue #	Diameter	Length
395451	.018"	40cm
395200	.018"	80cm
395513	.021"	40cm
395206	.025"	40cm

PTFE Coated Stainless Steel (Box of 10)

Catalogue #	Diameter	Length
395225	.025"	150cm
395687	.035"	150cm
395231	.035"	150cm

Exchange Wires - PTFE Coated (Box of 10)

Catalogue #	Diameter	Length
395367	.025"	260cm
395369	.035"	260cm
395370	.038"	260cm


Guidewires - Exchange Wires

Box of 10

Catalogue #	Diameter	Length	Tip Curve Radius
395585	.025"	260cm	3mm
395722	.035"	200cm	3mm
395794	.035"	260cm	3mm

Lunderquist Stainless Steel Guidewires

The Lunderquist Stainless Steel Guidewire is an extremely stiff guidewire used to provide extra support during non-vascular procedures. Flexible tip for access to delicate or difficult areas; Lunderquist Guidewire for firm support; available in a variety of lengths and tip configurations.



Lunderquist Stainless Steel Guidewire - Straight (Box of 10)

Catalogue #	Diameter	Length
105128082	.028"	82cm
105135082	.035"	82cm

Lunderquist Stainless Steel Guidewire - 3mm J Curve (Box of 10)


Catalogue #	Diameter	Length
105235082	.035"	82cm
105235120	.035"	120cm

POINTER™ Nitinol Guidewire

Flexibility and stability, strength and malleability all in one wire – available in size .018". Excellent visibility through highly radiopaque coil tip; smooth transition from nitinol core to hydrophilic coated radiopaque tip provides excellent manoeuvrability; hydrophilic coating on tip for less traumatic path finding/ advancement; PTFE-coated nitinol shaft for minimal friction during introduction and manipulation; Nitinol shaft construction provides 1:1 torque control and kink-resistance for negotiating through tortuous anatomy.

Box of 5

Catalogue #	Diameter	Length
180418210	.018"	210cm
180418260	.018"	260cm
180418300	.018"	300cm



Guidewires - Specialty

Stainless Steel (Box of 10)

Catalogue #	Diameter	Length	Tip Curve Radius
395573	.025"	40cm	3mm
395271	.035"	40cm	3mm
395654	.035"	50cm	3mm

Stainless Steel Double-ended (Box of 10)

Catalogue #	Diameter	Length	Tip Curve Radius
395539	.032"	65cm	Str & 3mm
395680	.035"	45cm	Str & 3mm
395899	.035"	65cm	Str & 3mm

Bentson Wires - PTFE Coated (Box of 10)

Catalogue #	Diameter	Length
395170	.035"	180cm
395930E	.035"	150cm

Guidewires - Movable Core

Movable Core PTFE Coated (Box of 10)

Catalogue #	Diameter	Length	Tip Curve Radius
395312	.035"	150cm	3mm

.018" Access Wires

Key Product Attributes:

- Smooth transition from rigid core to radiopaque coil tip provides excellent manoeuvrability
- Nitinol offers kink-resistance and torque control
- Stainless steel for maximum backup
- All guidewires have a flexible radiopaque tip

Box of 5

Catalogue #	Diameter	Length	Core Material
107018080	.018"	80cm	Stainless Steel
172618080	.018"	80cm	Nitinol



DUOMED

Duomed Swiss AG

Grenzstrasse 5a, 6214 Schenkon

Tel. +41 (0)41 510 07 00

www.duomed.com