



DUOMED

Duomed Swiss AG
Grenzstrasse 5a, 6214 Schenkon
Tel. +41 (0)41 510 07 00
www.duomed.com

Sets & Arterial Needles


- **V-Stick™ Vascular Access Sets**
- **Axcess™ Percutaneous Introducer Kits**
- **OneStic™ Micro-Access Introducer Kits**
- **Arterial Needles**
- **Guidewire Introducer Needles**
- **Introducer Kits for Venous Applications**



V-Stick™ Vascular Access Sets

V-Stick Vascular Access Sets assist in gaining vascular access for the placement of .035" and .038" guidewires into the vascular system using small needle access to reduce puncture site and vessel trauma.

V-Stick Vascular Access Sets include an UltraSharp 21G needle featuring a silicone-coated tip and a coaxial introducer, available in standard or stiff, designed to enable smooth transition.



V-Stick Vascular Access Sets w/Standard Coaxial Introducer (Box of 10)

Catalogue #	Sheath	Needle			Dilator		Guidewire		
	Size	Gauge	Length	Echogenic	Stiffness	Locking Nut Colour	Diameter	Length	Material
391504300E	4F	21G	7cm	Yes	Standard	Red	.018"	40cm	NT
391504400E	4F	21G	7cm	No	Standard	Red	.018"	40cm	SS
391504500E	4F	21G	7cm	No	Standard	Red	.018"	40cm	SS
391505000E	5F	21G	7cm	Yes	Standard	Red	.018"	40cm	SS
391505300E	5F	21G	7cm	Yes	Standard	Red	.018"	40cm	NT
391505400E	5F	21G	7cm	No	Standard	Red	.018"	40cm	NT
391505500E	5F	21G	7cm	No	Standard	Red	.018"	40cm	SS

SS = Stainless Steel


NT = Nitinol

V-Stick Vascular Access Sets w/Stiff Coaxial Introducer (Box of 10)

Catalogue #	Sheath	Needle			Dilator		Guidewire		
	Size	Gauge	Length	Echogenic	Stiffness	Locking Nut Colour	Diameter	Length	Material
391604000E	4F	21G	7cm	Yes	Stiff	White	.018"	40cm	SS
391604300E	4F	21G	7cm	Yes	Stiff	White	.018"	40cm	NT
391604400E	4F	21G	7cm	No	Stiff	White	.018"	40cm	NT
391604500E	4F	21G	7cm	No	Stiff	White	.018"	40cm	SS
391605000E	5F	21G	7cm	Yes	Stiff	White	.018"	40cm	SS
391605300E	5F	21G	7cm	Yes	Stiff	White	.018"	40cm	NT
391605400E	5F	21G	7cm	No	Stiff	White	.018"	40cm	NT

Axcess™ Percutaneous Introducer Kits

Axcess Percutaneous Introducer Kits feature a valved introducer to reduce procedural blood loss and air entry. In addition to its ease of use, our Axcess Percutaneous Introducers are available in various lengths and sizes and the kit is available with or without a guidewire.



Axcess™ Percutaneous Introducer Kit w/Double-ended Wire (Box of 10)

Catalogue #	Sheath			Guidewire	
	Size	Length	Hub Colour	Diameter	Length
496004	4F	11cm	Red	.035"	45cm
496005	5F	11cm	Grey	.038"	45cm
496006	6F	11cm	Green	.038"	45cm
496007	7F	11cm	Orange	.038"	45cm
496008	8F	11cm	Blue	.038"	45cm
496009	9F	11cm	Black	.038"	45cm

OneStic™ Micro-Access Introducer Kits

OneStic Micro-Access Introducer Kits combine Argon's micro puncture set and valved introducer into one convenient kit. The OneStic Micro-Access Introducer Kit allows for time and procedural step reduction by including a 21G needle, a .018 guidewire, and various size and length sheath options.



OneStic Micro-Access Introducer Kit w/Stainless Steel Mandrel Wires (Box of 10)

Catalogue #	Sheath			Guidewire			Needle		
	Size	Length	Hub Colour	Diameter	Length	Core/Tip Material	Gauge	Length	Echogenic
497850	6F	5.75cm	Green	.018"	45cm	SS/SS	21G	5.1cm	No
497820	4F	11cm	Red	.018"	45cm	SS/SS	21G	7cm	No
497830	5F	11cm	Grey	.018"	45cm	SS/SS	21G	7cm	No
497840	6F	11cm	Green	.018"	45cm	SS/SS	21G	7cm	No
497861	7F	11cm	Orange	.018"	45cm	SS/SS	21G	7cm	No
497862	8F	11cm	Blue	.018"	45cm	SS/SS	21G	7cm	No


OneStic Micro-Access Introducer Kit w/Nitinol Wires with Palladium Tip (Box of 10)

Catalogue #	Sheath			Guidewire			Needle		
	Size	Length	Hub Colour	Diameter	Length	Core/Tip Material	Gauge	Length	Echogenic
497865	5F	5.75cm	Grey	.018"	40cm	NT/PD	21G	2.5cm	No
497841	6F	11cm	Green	.018"	40cm	NT/PD	21G	7cm	No

Introducer Kits for Venous Applications

Kit Contents:

- (1) Radiopaque Sheath
- (1) Vessel Dilator
- (1) Guidewire (Double-ended, 3mm J, and ST)




Tuohy-Borst Introducer Kit w/Integral Hemostasis Valve (Box of 10)

Catalogue #	Sheath			Guidewire			
	Size	Length	Hub Colour	Diameter	Length	Material	Tip Shade
496132	6F	11cm	Green	.035"	45cm	SS	Double-ended, 3mm J & Straight
496133	7F	11cm	Orange	.035"	45cm	SS	Double-ended, 3mm J & Straight
496134	8F	11cm	Blue	.035"	45cm	SS	Double-ended, 3mm J & Straight
496135	8.5F	11cm	White	.035"	45cm	SS	Double-ended, 3mm J & Straight
496136	9F	11cm	Black	.035"	45cm	SS	Double-ended, 3mm J & Straight

SS = Stainless Steel NT = Nitinol PD = Palladium

Introducer Sets & Arterial Needles



PCI Kit/2 - Introducer Kit w/Detachable Hemostasis Valve (Box of 10)

Catalogue #	Sheath		Guidewire	
	Size	Length	Hub Colour	Compatibility
496124	8F	11cm	Orange	up to .038"

Kit Contents: (1) Vessel Dilator and (1) Introducer Assembly

PCI Kit/5 - Introducer w/Detachable Hemostasis Valve (Box of 10)

Catalogue #	Sheath			Guidewire			
	Size	Length	Hub Colour	Diameter	Length	Material	Tip Shade
496183	7F	11cm	Orange	.035"	45cm	SS	Double-ended, 3mm J & Straight

Kit Contents:

- (1) Disposable Towel
- (1) Vessel Dilator
- (3) 4" x 4" Gauze Sponges
- (1) Guidewire (Double-ended, 3mm J, and ST)
- (1) Sheath with Detachable Hemostasis Valve

SS = Stainless Steel

Kit Contents:

- (1) Radiopaque Sheath
- (1) Vessel Dilator
- (1) Guidewire (Double-ended, 3mm J, and ST)

**One-Piece Introducer Kit w/Integral Hemostasis Valve (Box of 10)**

Catalogue #	Sheath			Guidewire			
	Size	Length	Hub Colour	Diameter	Length	Material	Tip Shade
496150	4F	11cm	Red	.025"	45cm	SS	Double-ended, 3mm J & Straight
496151	5F	11cm	Grey	.035"	45cm	SS	Double-ended, 3mm J & Straight
496152	6F	11cm	Green	.035"	45cm	SS	Double-ended, 3mm J & Straight
496153	7F	11cm	Orange	.035"	45cm	SS	Double-ended, 3mm J & Straight
496154	8F	11cm	Blue	.035"	45cm	SS	Double-ended, 3mm J & Straight
496155	8.5F	11cm	White	.035"	45cm	SS	Double-ended, 3mm J & Straight
496156	9F	11cm	Black	.035"	45cm	SS	Double-ended, 3mm J & Straight

Arterial Needles

AMC/3 - Arterial Needle w/Removable Needle Adapter (Box of 10)

Catalogue #	Gauge	Length	Guidewire Compatibility
195312	18G	7cm	up to .038"
195311	21G	2.5cm	up to .021"

AMC/4 - Clear Flashback Needle (Box of 10)

Catalogue #	Gauge	Length	Guidewire Compatibility
195317	18G	3.8cm	up to .038"
195318	18G	7cm	up to .038"
195319	19G	1.5cm	up to .028"
195322	21G	2.5cm	up to .021"
195323	21G	5.1cm	up to .021"

AMC/4 - Winged Arterial Needle (Box of 10)

Catalogue #	Gauge	Length	Guidewire Compatibility	Special Features
195352	18G	7cm	up to .038"	N/A

Percutaneous Entry Needles (Box of 10)

Catalogue #	Description	Gauge	Length	Guidewire Compatibility
195326	Percutaneous Entry Needle(Two-Part Thin Wall Tapered Cannula and Inner Arterial Beveled Stylet)	18G	7.3cm	up to .038"

Arterial Extensions (Box of 25)

Catalogue #	Description
041107002A	Arterial Extensions 7Ft. PVC Extension with Stopcock and Male Luer Lock Adaptor



Guidewire Introducer Needles

Guidewire Introducer Needles are designed for anterior, single wall, arterial percutaneous punctures. The tapered luer hub allows for easy guidewire insertion.

Box of 25

Catalogue #	Gauge	Length	Guidewire Compatibility	Special Features
GW11802	18G	7cm	up to .038"	
GW11902	19G	7cm	up to .038"	
GW11802WW	18G	7cm	up to .038"	Clear window wall hub
BGWI1802	18G	7cm	up to .038"	Winged base plate
BGWI1902	19G	7cm	up to .035"	Winged base plate



Seldinger Needles

Seldinger Needles are designed for easy, atraumatic, percutaneous arterial punctures.

Box of 25

Catalogue #	Gauge	Length	Guidewire Compatibility
SN1802	18G	7cm	up to .038"
SN1902	19G	7cm	up to .038"



Duomed Swiss AG
 Grenzstrasse 5a, 6214 Schenkon
 Tel. +41 (0)41 510 07 00
www.duomed.com