



DUOMED

Duomed Swiss AG

Grenzstrasse 5a, 6214 Schenkon

Tel. +41 (0)41 510 07 00

www.duomed.com


Mammography

- **Repositionable Breast Localization Needles (BLN)**
- **Non-Repositionable Breast Localization Needles (BLN)**
- **Needle Guides**
- **Galactography Kits**



Homer BLN

The Homer BLN contains a strong, flexible J-curve wire that may be retracted and repositioned as often as necessary to achieve accurate placement.




Box of 10

Catalogue #	Gauge	Length	Wire
231030G	20G	3cm	3cm
231050G	20G	5cm	11cm
231075G	20G	7.5cm	13.5cm
231100G	20G	10cm	16.1cm
231125G	20G	12.5cm	18.5cm

Hawkins™ I BLN

The Hawkins I BLN is repositionable to ensure accurate placement and is available in two wire choices:

- FlexStrand™, a 19-strand cable for maximum flexibility with minimal risk of transection
- A traditional hardwire, which minimizes flex during placement




Box of 10

Catalogue #	Gauge	Length	Wire
251050	20G	5cm	7.8cm
251075	20G	7.5cm	10.3cm
251100	20G	10cm	12.8cm
251125	20G	12.5cm	15.3cm

Hawkins™ II BLN

The Hawkins II BLN has a retention barb designed to provide strong anchoring strength in the breast.

- Wire is marked in three places: I Indicates wire at distal tip of needle; II Indicates side barb deployment; and III indicates distal barb deployment
- FlexStrand, a 19-strand cable for maximum flexibility with minimal risk of transection
- A traditional hardwire, which minimizes flex during placement



Hawkins II Breast Localization Needles w/FlexStrand (Box of 10)


Catalogue #	Gauge	Length	Wire
252050	20G	5cm	14.4cm
252075	20G	7.5cm	19.4cm
252100	20G	10cm	24.4cm
252125	20G	12.5cm	29.4cm

Hawkins II Breast Localization Needles w/Hardwire (Box of 10)

Catalogue #	Gauge	Length	Wire
242050	20G	5cm	14.4cm
242075	20G	7.5cm	19.4cm
242100	20G	10cm	24.4cm
242125	20G	12.5cm	29.4cm

Hawkins™ III BLN

The Hawkins III BLN has a barb design to provide strong anchoring strength in the breast. Wire is marked in two places: I indicates wire at distal tip of needle; II indicates distal barb deployment.



Hawkins III Breast Localization Needles w/FlexStrand™ (Box of 10)


Catalogue #	Gauge	Length	Wire
253030	20G	3cm	10.4cm
253050	20G	5cm	14.4cm
253075	20G	7.5cm	19.4cm
253100	20G	10cm	24.4cm
253125	20G	12.5cm	29.4cm

Hawkins III Breast Localization Needles w/Hardwire (Box of 10)

Catalogue #	Gauge	Length	w/Echogenic Tip	Wire
243030	20G	3cm	244030	10.4cm
243050	20G	5cm	244050	14.4cm
243075	20G	7.5cm	244075	19.4cm
243100	20G	10cm	244100	24.4cm
243125	20G	12.5cm	244125	29.4cm

Accura™ BLN (Kopans Style)

Accura BLN is a traditional, springhook hardwire with stiffener design. Wire is marked in two places: I indicates wire at distal tip of needle; II indicates distal barb deployment.




Accura Breast Localization Needles (Kopans Style) (Box of 10)

Catalogue #	Gauge	Length	Wire
BLN2003	20G	3cm	10.4cm
BLN2005	20G	5cm	14.4cm
BLN20075	20G	7.5cm	19.4cm
BLN2010	20G	10cm	24.4cm
BLN2015	20G	15cm	28.9cm
BLN2105	21G	5cm	14.4cm
BLN21075	21G	7.5cm	19.4cm
BLN2110	21G	10cm	24.4cm

Accura™ II BLN (Sadowsky Style)

Accura II BLN is designed for a preoperative marking of nonpalpable breast lesions. Wire is marked in two places: I indicates wire at distal tip of needle; II indicates distal barb deployment.




Accura II Breast Localization Needles (Sadowsky Style) (Box of 10)

Catalogue #	Gauge	Length	Stiffening Cannula	Wire
266050	20G	5cm	23G	14.4cm
266075	20G	7.5cm	23G	19.4cm
266100	20G	10cm	23G	24.4cm

D Wire BLN

The D Wire BLN is designed to provide strong anchoring strength in the breast.




D Wire Breast Localization Needles (Box of 10)

Catalogue #	Gauge	Length	Wire
BLND2005	20G	5cm	17.2cm
BLND20075	20G	7.5cm	18.3cm
BLND2010	20G	10cm	21.2cm

Needle Guides

Needle guides are designed specifically for use in conjunction with the Lorad Stereo Loc™ II Prone Stereotactic System for core biopsy of the breast.



Needle Guides for Lorad Stereo Loc II (Box of 20)

Catalogue #	Type	Compatible Device	Guide Quantity
LNGC14GX	Closed type	14G	1
LNGC1619PX	Closed type	16G to 19G	2
LNGC2022YX	Closed type	20G to 22G	2
LNGO1619PX	Open type	16G to 19G	2
LNGO2022YX	Open type	20G to 22G	2

Galactography Kits

Galactography kits are designed to perform safe and effective contrast media-enhanced radiography of the mammary ducts.



Galactography Kits (Box of 10)

Catalogue #	Gauge	Cannula
DGK002SX	30G	Straight
DGK002X	30G	Curved